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
Figure 1. Greater than 95% CHO Cell Transfection Efficiency using MaxCyte Transient Transfection. CHOs cells were transfected by static, small scale electroporation (EP) in an OC-400 processing assembly with a plasmid encoding green fluorescent protein (2µg DNA/1E6 cells). GFP expression and viability were measured by flow cytometry (FACS) 24 hrs post EP.

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
Rapid, Gram Scale Antibody Production

	Culture Volume	EP Volume	# of Cells	[IgG]	Total IgG Produced
Small Scale	20mL	0.4mL	8E7	1.40 g/L	28mg
Large Scale	2.8L	50mL	1E10	1.22 g/L	3.42 g

Figure 2: High Titer Antibody Production with Seamless Scalability.

8E7 or 1E10 CHO-S cells were transfected with an antibody expression plasmid (2 μ g DNA per 1E6 cells) via small scale (static), or large scale (flow) electroporation on the MaxCyte STX. Cells were plated at 6E6 cells/mL post electroporation. 1mM sodium butyrate was added to cultures and the temperature lowered to 32° C 24 hours post electroporation. Cultures were fed daily with a media optimized for antibody production. Total IgG was measured using ELISA on various days post transfection.

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
Figure 3: Similarity of Antibodies Produced from Stably and Transiently Transfected CHO Cells. CHO suspension cells were transfected by static electroporation with a plasmid encoding IgG heavy and light chain sequences. Secreted antibodies from electroporated cells (S) and a reference stable cell line (R) were examined by reducing and non-reducing SDS PAGE and Isoelectric Focusing.
Data courtesy of NovImmune.

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Titers of 1.2 grams/L Upon Culture Feeding Optimization


Figure 4: Substantially Higher Antibody Titers with Optimized Media Composition and Feeding Schedule.

8E7 CHO-S cells were resuspended at 2E8 cells/mL and transiently transfected with a human antibody expression plasmid (2µg DNA/1E6 cells) via static electroporation. Post electroporation cells were seeded at 5.6E6 cells/mL in Media 1 or Media 2. Cells were fed on two different schedules. Antibody titers were determined on days 3-14 post electroporation using an ELISA.

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
Figure 5: Post Transfection Inoculation at Higher Cell Densities Increases Antibody Titers.

1E10 CHO-S cells were transiently transfected with an antibody expression plasmid (1 μ g DNA/1E6 cells) via flow electroporation. Post electroporation, cells were inoculated into two shake flasks at densities of 6E6 or 1E7 cells per mL.

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Antibody Titers Maintained During Process Scale Up


Figure 6: Titers Maintained Upon Transfection and Culture Vessel Scale Up.

1E10 or 2.4E8 CHO-S cells were transfected with an antibody expression plasmid (1 μ g DNA/1E6 cells) using flow or static electroporation, respectively. Cells transfected via static EP were inoculated into 50mL of media in a 250mL shake flask. Cells transfected via flow EP were split and inoculated into 1L of media in a 3L shake flask or into 50mL of media in a 250mL shake flask. Starting cell densities were identical in all three flasks.

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Flow Electroporation for Stable Cell Line Generation


Figure 7: Rapid Stable Cell Development Using MaxCyte Electroporation.

A humanized mAb DNA plasmid was transfected into CHO-S cells using the CHO2 electroporation protocol on the MaxCyte STX. Limiting dilution cloning was carried out from the stable pool in G418 selection media. Cell lines were generated at week 6 and scaled up for production and accession cell bank. Productivity of the selected clone was > 3g/L as a fed batch.

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1. Strategies for Rapid Production of Therapeutic Proteins in Mammalian Cells. *BioProcess Int.* 10(4), 2012, p32-48.
2. Recombinant Protein Production by Large-scale Transient Gene Expression in Mammalian Cells: State of the Art and Future Perspectives. *Biotechnol. Lett.* 29(5), 2007, p677-684.
3. High-level Protein Expression in Scalable CHO Transient Transfection. *Biotechnol. Bioeng.* 103(3), 2009, p542-551.
4. Transient Transfection Factors for High-Level Recombinant Protein Production in Suspension Cultured Mammalian Cells. *Mol. Biotechnol.* 39(2), 2008, p141-153.
5. A High-Yielding, CHO-K1-Based Transient Transfection System. *BioProcess Int.* 11(1), 2013, p28-35.
6. Scalable Transient Gene Expression in Chinese Hamster Ovary Cells in Instrumented and Non-instrumented Cultivation Systems. *Biotechnol. Lett.* 29(5), 2007, p703-711.
7. Control of Culture Environment for Improved Polyethylenimine-Mediated Transient Production of Recombinant Monoclonal Antibodies by CHO Cells. *Biotechnol. Progr.* 22(3), 2006, p753-762.
8. Design of Experiment in CHO and HEK Transient Transfection Condition Optimization. *Prot. Expr. and Purif.* 78(1), 2011, p61-68.
9. Mild Hypothermia Improves Transient Gene Expression Yields Several Fold in Chinese Hamster Ovary Cells. *Biotechnol. Prog.* 24(2), 2008, p458-465.
10. Rapid Protein Production Using CHO Stable Transfection Pools. *Biotechnol. Progr.* 26(5), 2010, p1431-1437.